



Contents

1 Overview	3
2 Uniform Resource Identifiers.....	3
3 Device Management Interfaces	3
3.1 Update an existing patient	3
3.2 Remove device from patient	7

1 Overview

This AirView Exchange (AVX) technical specification only covers Device Management APIs exposed by AVX, which are intended for use by integration partners.

Before you can use the APIs outlined in this document, you must first acquire the patient's EasyCare Number (ECN). Refer to the *AVX API Technical Specification for Patient Management* for instructions on how to get a patient's ECN.

The information in this technical specification supplements the *AVX API Technical Specification for Patient Management*. It aligns with Release 4.6 of AirView for the Americas and the Asia-Pacific region.

2 Uniform Resource Identifiers

The following table depicts the hierarchical structured design of AVX information. It also summarizes the HTTP calls that are available at each level.

URIs are specified in a manner that is consistent with RFC 6570 URI Template, and are case sensitive.

URI	Method	Description
https://<host>/v3/patients/{ecn}	PUT	Used to update a patient record. See "3.1 Update an existing patient" on page 3.
https://<host>/v1/patients/{ecn}/device/{serialNumber}	DELETE	This removes a device with the provided serial number from the patient with the provided ECN. See "3.2 Remove device from patient" on page 7.

3 Device Management Interfaces

The following interface methods allow patients to be edited (to reassign a device), and remove devices from patients.

3.1 Update an existing patient

This API call is used to update an existing patient record in AirView. The patient must be accessible by the IA.

Note: This API is also referenced in section 7.3.2 of the *AVX API Technical Specification for Patient Management*. It is required to update a patient's device serial number.

This version supports S9™, Air10™, Stellar™ and Astral™ devices only. It is also valid for a patient to not have a device.

Property	Value
Method	PUT
URL	https://<host>/v3/patients/{ecn}
	where: <i>ecn</i> is the unique identifier of the patient being updated.
PUT Content	A "Patient" JSON object – See table below

The table on the following page identifies the mandatory and optional fields contained within a "patient" JSON object for a patient create. If a field is not included, it is preserved.

JSON Name	Type	Req'd	Rules/Format	Description
firstName	string	No	1 to 50 characters	Patient's first name
lastName	string	No	1 to 50 characters	Patient's last name
dob	string	No	YYYY-MM-DD	Patient's date of birth
location	string	No	Valid location name in organization	The location name of the patient's assigned location. See "Retrieve an organization's locations" in the AVX API Technical Specification for Patient Management.
patientID	string	No	1 to 20 characters	Patient identifier
clinicalUser	string	No	Valid clinical user ID within organization	The user ID of the patient's assigned clinical user. See "Retrieve clinical users who have access to a given location in an organization" in the AVX API Technical Specification for Patient Management.
setupDate	string	No	YYYY-MM-DD	The date from which a patient's wireless monitoring schedule begins and from which compliance calculation/s begin/s with format YYYY-MM-DD
physician	string	No	Valid physician user ID	The physician's user ID is included as an associated physician to the created patient. See "Get physician list" in the AVX API Technical Specification for Patient Management.
payorId	string	No	Valid payor ID	Payor ID to be associated with the patient. If it is not included, the organization's default payor is assigned. See "Get payor list" in the AVX API Technical Specification for Patient Management.
prescription ()	object	No	n/a	Patient prescription information
module ()	object	No ¹	n/a	Module object
serialNo	string	No	Valid module serial number allocated to org and not assigned.	Patient's communications module serial number
device ()	object	No	n/a	Device object
serialNo	string	No ²	Valid device serial number not assigned.	Patient's prescribed device serial number
deviceType	string	No ²	See "S9 devices" in the AVX API Technical Specification for Patient Management.	Patient's prescribed device type. Note: This is applicable for S9 devices only.
deviceNumber	string	No ³	3-digit number	Device number is required for Air10, Stellar, AirCurve™ 10 ST-A and Astral devices.

¹ module() is applicable only for an S9 device, and must not be sent for other devices.

² deviceType is required if a device is supplied and it is an S9. deviceType for Air10, Stellar, AirCurve 10 ST-A and Astral devices are known by AirView and as such must not be sent.

³ deviceNumber is required for Air10, Stellar, AirCurve 10 ST-A and Astral devices.

Example (for S9):

```
{
  "patient": {
    "dob": "1970-01-01",
    "firstName": "General",
    "lastName": "Hauck",
    "location": "Unimedis AL",
    "clinicalUser": "28a2898d-4e8e-45df-ae08-526337b74920",
    "setupDate": "1990-01-01",
    "physician": "c501e66d-bcd0-4d99-b745-0970e65d55dd",
    "prescription": {
      "module": {
        "serialNo": "09876543210"
      },
      "device": {
        "serialNo": "01234567890",
        "deviceType": "026001"
      }
    }
  }
}
```

Note: Identities and data are fictional.

Example (for Air10, Stellar and Astral devices):

```
{
  "patient": {
    "dob": "1972-01-01",
    "firstName": "General",
    "lastName": "Williamson",
    "location": "Unimedis AP",
    "patientID": "UTY766555",
    "clinicalUser": "c501e66d-bcd0-4d99-b745-0970e65d55dd",
    "setupDate": "1999-01-01",
    "physician": "28a2898d-4e8e-45df-ae08-526337b74920",
    "prescription": {
      "device": {
        "serialNo": "12547854521",
        "deviceNumber": "132"
      }
    }
  }
}
```

Note: Identities and data are fictional.

Responses:

Success:

Property	Value
HTTP Response Status Code	200
Content Type	application/json
Body	JSON "Response" Object – See table below

The following fields are included in the JSON "Response" Object.

Field Name (JSON Name)	Type	Description
responseCode	string	1002000
responseDesc	string	Unspecified – Debugging purposes only
ecn	string	Unique identifier for the updated patient

Example:

```
{
  "Response": {
    "responseCode": "1002000",
    "responseDesc": "Patient Data Updated",
    "ecn": "44187887-a48d-4f4e-9393-e4207045a2c4"
  }
}
```

Note: Data is fictional.

Errors:

Property	Value
HTTP Response Status Code	Function of error – See table below
Content Type	application/json
Body	JSON "Response" Object – See table below

Field Name (JSON Name)	Type	Description
responseCode	string	Function of error – See table below
responseDesc	string	Unspecified – Debugging purposes only
fieldName	string	Any applicable field name – Debugging purposes only

Error Description	Response Type	responseCode
Requesting user does not have access to the patient	HTTP 401	1002001
Mandatory Field Missing	HTTP 400	1002002
JSON object does not match schema	HTTP 400	1002010
Invalid Parameter General	HTTP 400	1001003
Device serial number already assigned to another patient	HTTP 403	1002101
Module serial number already assigned to another patient	HTTP 403	1002102
Module serial number not assigned to your organization	HTTP 403	1002103
Device type included without a device serial number	HTTP 400	1002100
Device serial number included without device type	HTTP 400	1002100
Incorrect serial number/device number combination	HTTP 403	1002105

3.2 Remove device from patient

This API call can be used to remove a device from a particular patient. Note: This version of the API cannot be used to remove an S9 module.

Property	Value
Method	DELETE
URL	https://<host>/v1/patients/{ecn}/device/{serialNumber}
	where: <i>ecn</i> is the patient unique identifier. <i>serialNumber</i> is the device serial number to be removed from the patient file.

Responses:

Success:

An HTTP 200 response with the following fields:

Property	Value
HTTP Response Status Code	200
Content Type	application/json
Body	JSON "Response" Object – See table below

Errors:

Property	Value
HTTP Response Status Code	Function of error – See table below
Content Type	application/json
Body	JSON "Response" Object – See table below

Field Name (JSON Name)	Type	Description
responseCode	string	Function of error – See table below
responseDesc	string	Unspecified – Debugging purposes only

Description	Response Type	Error Code
Requesting user does not have access to the patient	HTTP401	1002001
Specified serial number does not belong to patient	HTTP404	1013001



ResMed Corp

9001 Spectrum Center Boulevard San Diego CA 92123 USA

See www.resmed.com for other ResMed locations worldwide.

For patent and other intellectual property information, see www.resmed.com/ip.

TransUnion is a registered trademark of Trans Union LLC. AcuCare, Air10, AirCurve, AirFit, AirMini, AirSense, AirView, Astral, Mirage, Mirage Activa, Mirage Liberty, Mirage Quattro, Mirage Swift, PIXI, Quattro, S9, Stellar, Swift, and Ultra Mirage are trademarks and/or registered trademarks of the ResMed family of companies. 2281354/3 2021-06

ResMed.com