



Contents

1 Overview	3
2 Authentication	3
3 Integrator types	3
4 Match patients	3
4.1 Workflow overview	3
4.2 Workflow tips	3
4.3 Required APIs	3
4.4 Instructions.....	4
5 Create or update patients	4
5.1 Workflow overview	4
5.2 Workflow tips	4
5.3 Required APIs	4
5.4 Instructions.....	4
6 Get patient details	5
6.1 Workflow overview	5
6.2 Workflow tips	5
6.3 Required APIs	5
6.4 Instructions.....	5

1 Overview

AirView Exchange (AVX) (formerly known as ResMed EasyCare Exchange/ECX2) is a complete integration Application Programming Interface (API) suite designed for two-way communication from certified applications into the ResMed AirView system.

This document covers best practices for the optimisation of patient management APIs so you can design your integration with efficiency in mind. We highly recommend that you review each workflow to see what best practices you can apply to your integration setup.

2 Authentication

Before you continue, ensure your integrating application is correctly authenticated with AVX. You must include the authentication credentials for each call to AVX. For more information, refer to section 3 *Architecture overview* in the *AVX API Technical Specification for Patient Management*.

3 Integrator types

AVX supports multiple integrator types. Most accounts are set up with one of the following integrator types:

Organisation Specific Integrator: You have access to all patients in a HME's AirView organisation. A HME administrator must provide access to your integrator account in AirView.

Patient Level Integrator: You only have access to patients that the HME has granted you permission to access. A HME clinical user can do this while they set up a new patient in AirView or edit an existing patient in their organisation.

4 Match patients

4.1 Workflow overview

The instructions below outline how you can set up a patient matching framework to associate patients in your system with patients in AirView.

Follow this workflow when you set up your integrating application. You need:

- a patient matching framework to associate your system's patient ID with AirView's patient ID228
- to reconcile your patient list in AirView if you did not create patients with your integrating application through AVX
- to have your integration recognise which ECNs were already processed.

Note: If you have pre-existing patients in AirView, but you want to manage the creation of AirView patients going forward, you need to complete this workflow.

4.2 Workflow tips

Your application must use **Get Accessible Organisations** regardless of what kind of integration it is. You can have multiple organisations assigned to each integrator type.

4.3 Required APIs

- **Get Accessible Organisations**
(GET <host>/v1/orgs)
- **Get Patient List**
(GET <host>/v1/orgs/{orgID}/patients)
- **Get Patient Details**
(GET <host>/v1/patients/{ecn})

4.4 Instructions

1. Call the **Get Accessible Organisations** API to get a list of organisations you have access to (this may only be one organisation ID).
2. With the organisations ID you retrieved from step 1, call the **Get Patient List** API to get a list of ECNs in your organisation.
Note: We recommend that you save this list in your patient database. If you are set up as a patient level integrator, you will only receive ECNs that you have access to (see "3 Integrator types" on page 3).
3. Use the **Get Patient Details** API for each retrieved ECN to get the patient's first name, last name and date of birth.
4. Use the patient's first name, last name and date of birth to find a match in your system.
5. Once you find a match, associate the patient's ECN with your system's unique ID.

5 Create or update patients

5.1 Workflow overview

The instructions below outline how you can create or update a patient in your system with AVX.

Follow this workflow if your integrating application creates patients directly in AirView.

5.2 Workflow tips

- The Patient ID field in AirView, which is different from the ECN, is an optional field that you can add when you create or update a patient record in AirView.
- In special-use cases, you should cache the list of organisations and locations. For example, when you configure an available subset per user or relatively static organisation structure. If you cache, your application must refresh its cache when the organisation structure changes because of an administrative action.
- If you want to update the patient's demographic information, patientId, setupDate or prescription details, you should not call to get the organisation and location. You can use the saved ECNs from your system to update the patient's demographic information. Only use these calls to update the patient's location.

5.3 Required APIs

- **Get Accessible Organisations**
(GET <host>/v1/orgs)
- **Get Locations for Organisation**
(GET <host>/v1/orgs/{orgID}/locations[?index={index}&count={count}])
- **Create Patient** (Create a new patient at a specific location for a clinical user)
(POST <host>/v3/orgs/{orgID}/patients)

5.4 Instructions

1. Call the **Get Accessible Organisations** API to get a list of organisations you have access to (this may only be one organisation ID).
2. For each organisation, call the **Get Locations for Organisation** API.
Note: We recommend that you save the location list in your integrating application and check for updates every two weeks.
3. Call the **Create Patient** API to create a patient and assign them to a location.
Note: Alternatively, you can update the patient at this step.

When you complete this procedure, the API returns an ECN for the newly created patient. We recommend that you store this ECN to access other AVX APIs in the future.

Note: When the API displays the locations, call the **Get Locations for Organisation** API only when you select a specific organisation.

6 Get patient details

6.1 Workflow overview

The instructions below outline how you can retrieve a patient's details from AirView.

Follow this workflow to retrieve a patient's demographics and basic device information.

Note: You can only use this workflow if you store patient ECNs in your system. We recommend that you only use this API once every one to two weeks.

6.2 Workflow tips

- You can always find a patient's ECN at the end of the AirView URL.
For example, `https://<host>/patients/{ecn}`
- The **Get Patient List** API does not return ECNs in any order.
- If you want to keep your system's patient demographic details up to date with changes in AirView, we recommend you call the **Get Patient Details** API once every two weeks or monthly as these details are not updated often.

6.3 Required APIs

- **Get Accessible Organisations**
(GET <host>/v1/orgs)
- **Get Patient List**
(GET <host>/v1/orgs/{orgID}/patients)
- **Get Patient Details**
(GET <host>/v1/patients/{ecn})

6.4 Instructions

Use the stored ECNs from your database to call the **Get Patient Details** API.



ResMed Corp

9001 Spectrum Center Boulevard San Diego CA 92123 USA

See www.resmed.com for other ResMed locations worldwide.

For patent and other intellectual property information, see www.resmed.com/ip.

AirView is a trademark and/or registered trademark of the ResMed family of companies. 2281471/3 2021-01

ResMed.com