



Contents

1	Overview.....	3
2	Uniform Resource Identifiers.....	3
3	Patient note interfaces.....	3
3.1	Create a new patient note	3
3.2	Get a patient note.....	6
3.3	Get patient notes by date range	8

1 Overview

This AirView Exchange (AVX) technical specification only covers Patient Notes APIs exposed by AVX, which are intended for use by integration partners.

Before you can use the APIs outlined in this document, you must first acquire the patient's EasyCare Number (ECN). Refer to the *AVX API Technical Specification for Patient Management* for instructions on how to get a patient's ECN.

The information in this technical specification supplements the *AVX API Technical Specification for Patient Management*. It aligns with Release 4.11 of AirView for the Americas and the Asia-Pacific region.

2 Uniform Resource Identifiers

The following table depicts the hierarchical structured design of AVX information. It also summarizes the HTTP calls that are available at each level.

URIs are specified in a manner that is consistent with RFC 6570 URI Template, and are case sensitive.

URI	Method	Description
https://<host>/v1/patients/{ecn}/notes	POST	Used to create a new patient note for the specified ECN. See "3.1 Create a new patient note" on page 3.
https://<host>/v1/patients/{ecn}/notes/{notedId}	GET	Returns a patient note for the requested ECN and note identifier. See "3.2 Get a patient note" on page 6.
https://<host>/v1/patients/{ecn}/notes?startDate={startDate}&endDate={endDate}	GET	Returns a list of patient notes for the requested ECN for a defined date range. See "3.3 Get patient notes by date range" on page 8.

3 Patient note interfaces

The following interface methods allow new patient notes to be created and existing patient notes to be retrieved.

3.1 Create a new patient note

This API call can be used to create a new patient note in AirView for a patient accessible by the Integrating Application (IA). An AirView user may be specified as the creator of the note. If an AirView user is not specified, then the name and organization of an IA user must be specified as the creator of the note.

Property	Value
Method	POST
URL	https://<host>/v1/patients/{ecn}/notes where: <i>ecn</i> is the unique identifier of the patient. example: https://<host>/v1/patients/12345678-11AA-22BB-33CC-1234567890AA/notes
POST Content	JSON "createPatientNote" object – see table below.

The following table identifies the contents of the JSON “createPatientNote” object. It is important to verify the proper case and spelling of optional field names. Misspelt optional field names will be ignored.

JSON Name	Type	Required	Rules/Format	Description
body	string	Yes	1 to 400 characters	The text content of the note. It can't be empty or contain only whitespace characters.
createdDate	string	Yes	ISO 8601 combined date and time extended format	The date and time the note will be listed as created on. It must not be later than the date and time the request was received by the server.
createdUserId	string	No*	Valid AirView user ID	The unique identifier of the AirView user requesting to create the note.
externalUserOrganization	string	No*	1 to 50 characters	The organization of the IA user requesting to create the note.
externalUserFirstName	string	No*	1 to 50 characters	The first name of the IA user requesting to create the note.
externalUserLastName	string	No*	1 to 50 characters	The last name of the IA user requesting to create the note.

* If createdUserId is provided, then externalUserFirstName, externalUserLastName and externalUserOrganization are optional, otherwise they are required.

Example (without createdUserId):

```
{
  "createPatientNote": {
    "body": "The patient called with a support question",
    "createdDate": "2018-11-19T10:10:04-04:00",
    "externalUserOrganization": "Large HME",
    "externalUserFirstName": "Jane",
    "externalUserLastName": "Doe"
  }
}
```

Note: Data is fictional.

Example (with createdUserId):

```
{
  "createPatientNote": {
    "body": "The patient called with a support question",
    "createdDate": "2018-11-19T14:10:04Z",
    "createdUserId": "1234"
  }
}
```

Note: Data is fictional.

Responses:

Success:

Property	Value
HTTP Response Status Code	201
Content Type	application/json
Body	JSON “Response” object – see table below.

The following fields will be included in the JSON "Response" object:

Field Name (JSON Name)	Type	Description
responseCode	string	1002000
responseDesc	string	Unspecified - Debugging purposes only
patientNoteld ()	object	Contains the identifier of the newly added patient note.
id	string	The unique identifier of the newly added note.

Example:

```
{
  "Response": {
    "responseCode": "1002000",
    "patientNoteId": {
      "id": "12"
    }
  }
}
```

Note: Data is fictional.

Errors:

Property	Value
HTTP Response Status Code	See table below.
Content Type	application/json
Body	JSON "Response" object – see table below.

The following fields will be included in the JSON "Response" object:

Field Name (JSON Name)	Type	Description
responseCode	string	See table below.
responseDesc	string	Unspecified - Debugging purposes only
fieldName	string	Any applicable field name – Debugging purposes only

The following response code values may be returned:

Error Description	HTTP Response Status Code	responseCode
The request entity was not valid JSON.	HTTP 400	1002010
- The request entity was valid JSON but did not conform to the request schema defined above. - The createdDate provided was in the future. - The createdUserId was not found.	HTTP 400	1002001
- The client is not authorized to make this request. - The requested patient was not accessible or was not found.	HTTP 401	1002001
A note already exists for the requested patient and createdDate.	HTTP 422	1002001

Example:

```
{
  "Response": {
    "responseCode": "1002001"
  }
}
```

Note: Data is fictional.

3.2 Get a patient note

This API call can be used to retrieve a patient note for a patient accessible by the IA.

Property	Value
Method	GET
URL	https://<host>/v1/patients/{ecn}/notes/{noteld}
where: <i>ecn</i> is the unique identifier of the patient. <i>noteld</i> is the unique integer identifier of the note. The noteld can be obtained when a new note is created (see "3.1 Create a new patient note" on page 3) or when the list of patient notes is retrieved by date range (see "3.3 Get patient notes by date range" on page 8). example: https://<host>/v1/patients/12345678-11AA-22BB-33CC-1234567890AA/notes/12	

Responses:

Success:

Property	Value
HTTP Response Status Code:	200
Content Type:	application/json
Body:	JSON "Response" object – see table below.

The following fields will be included in the JSON "Response" object:

Field Name (JSON Name)	Type	Description
responseCode	string	1002000
responseDesc	string	Unspecified - Debugging purposes only
patientNote ()	object	Contains the requested patient note.
id	string	The unique identifier of the note.
body	string	The text content of the note.
createdDate	string	The date and time the note was created, in ISO 8601 combined date and time extended format.
createdUserId	string	The unique identifier of the AirView user who created the note. This is optionally present.

Field Name (JSON Name)	Type	Description
createdUserFirstName	string	The first name of the AirView user who created the note. This is optionally present; it is not present if createdUserId is not present or if the user has been removed from the system.
createdUserLastName	string	The last name of the AirView user who created the note. This is optionally present; it is not present if createdUserId is not present or if the user has been removed from the system.
externalUserOrganization	string	The organization of the IA user who created the note. This is optionally present.
externalUserFirstName	string	The first name of the IA user who created the note. This is optionally present.
externalUserLastName	string	The last name of the IA user who created the note. This is optionally present.

Example (without createdUserId):

```
{
  "Response": {
    "responseCode": "1002000",
    "patientNote": {
      "id": "12",
      "body": "The patient called with a support question",
      "createdDate": "2018-11-19T14:10:47Z",
      "externalUserOrganization": "Large HME",
      "externalUserFirstName": "Jane",
      "externalUserLastName": "Doe"
    }
  }
}
```

Note: Data is fictional.

Example (with createdUserId):

```
{
  "Response": {
    "responseCode": "1002000",
    "patientNote": {
      "id": "56",
      "body": "The patient called with a support question",
      "createdDate": "2018-11-19T14:10:47Z",
      "createdUserId": "1234",
      "createdUserFirstName": "John",
      "createdUserLastName": "Smith"
    }
  }
}
```

Note: Data is fictional.

Errors:

Property	Value
HTTP Response Status Code	See table below.
Content Type	application/json
Body	JSON "Response" object – see table below.

The following fields will be included in the JSON “Response” object:

Field Name (JSON Name)	Type	Description
responseCode	string	See table below.
responseDesc	string	Unspecified - Debugging purposes only
fieldName	string	Any applicable field name – Debugging purposes only

The following response code values may be returned:

Error Description	HTTP Response Status Code	responseCode
- The client is not authorized to make this request. - The requested patient was not accessible or was not found.	HTTP 401	1002001
The requested note was not accessible or was not found.	HTTP 404	1002001

Example:

```
{
  "Response": {
    "responseCode": "1002001"
  }
}
```

Note: Data is fictional.

3.3 Get patient notes by date range

This API call can be used to retrieve a list of patient notes for a patient accessible by the IA where the note `createdDate` falls within the requested date range. The list of notes is sorted in descending order by the `createdDate`, so the newer notes appear before the older notes.

Property	Value
Method	GET
URL	<code>https://<host>/v1/patients/{ecn}/notes?startDate={startDate}&endDate={endDate}</code> <i>where:</i> <i>ecn</i> is the unique identifier of the patient. <i>startDate</i> is the start date and time of the requested date range, in ISO 8601 combined date and time extended format. <i>endDate</i> is the end date and time of the requested date range, in ISO 8601 combined date and time extended format. <i>example:</i> <code>https://<host>/v1/patients/12345678-11AA-22BB-33CC-1234567890AA/notes?startDate=2018-01-01T00:00:00Z&endDate=2018-11-20T03:00:00Z</code> <code>https://<host>/v1/patients/12345678-11AA-22BB-33CC-1234567890AA/notes?startDate=2018-01-01T00:00:00Z&endDate=2018-11-20T14:00:00-04:00</code>

Responses:

Success:

Property	Value
HTTP Response Status Code	200
Content Type	application/json
Body	JSON "Response" object – see table below.

The following fields will be included in the JSON "Response" object:

Field Name (JSON Name)	Type	Description
responseCode	string	1002000
responseDesc	string	Unspecified - Debugging purposes only
patientNotes ([[]])	object array	Contains the requested patient notes.
id	string	The unique identifier of the note.
body	string	The text content of the note.
createdDate	string	The date and time the note was created, in ISO 8601 combined date and time extended format.
createdUserId	string	The unique identifier of the AirView user who created the note. This is optionally present.
createdUserFirstName	string	The first name of the AirView user who created the note. This is optionally present; it is not present if createdUserId is not present or if the user has been removed from the system.
createdUserLastName	string	The last name of the AirView user who created the note. This is optionally present; it is not present if createdUserId is not present or if the user has been removed from the system.
externalUserOrganization	string	The organization of the IA user who created the note. This is optionally present.
externalUserFirstName	string	The first name of the IA user who created the note. This is optionally present.
externalUserLastName	string	The last name of the IA user who created the note. This is optionally present.

Example (without createdUserId):

```
{
  "Response": {
    "responseCode": "1002000",
    "patientNotes": [
      {
        "id": "12",
        "body": "The patient called with a support question",
        "createdDate": "2018-11-19T22:10:03Z",
        "externalUserOrganization": "Large HME",
        "externalUserFirstName": "Jane",
        "externalUserLastName": "Doe"
      },
      {
        "id": "34",
        "body": "The patient called with a support question",
        "createdDate": "2018-11-19T14:10:47Z",
        "externalUserOrganization": "Large HME",
        "externalUserFirstName": "John",
        "externalUserLastName": "Doe"
      }
    ]
  }
}
```

Note: Data is fictional.

Example (with createdUserId):

```
{
  "Response": {
    "responseCode": "1002000",
    "patientNotes": [
      {
        "id": "56",
        "body": "The patient called with a support question",
        "createdDate": "2018-11-19T22:10:03Z",
        "createdUserId": "1234",
        "createdUserFirstName": "Jane",
        "createdUserLastName": "Smith"
      },
      {
        "id": "78",
        "body": "The patient called with a support question",
        "createdDate": "2018-11-19T14:10:47Z",
        "createdUserId": "5678",
        "createdUserFirstName": "John",
        "createdUserLastName": "Smith"
      }
    ]
  }
}
```

Note: Data is fictional.

Errors:

Property	Value
HTTP Response Status Code	See table below.
Content Type	application/json
Body	JSON "Response" object – see table below.

The following fields will be included in the JSON "Response" object:

Field Name (JSON Name)	Type	Description
responseCode	string	See table below.

Field Name (JSON Name)	Type	Description
responseDesc	string	Unspecified - Debugging purposes only
fieldName	string	Any applicable field name – Debugging purposes only

The following response code values may be returned:

Error Description	HTTP Response Status Code	responseCode
- The startDate or endDate parameters were missing or not in the expected ISO 8601 combined date and time extended format. - The endDate occurs before startDate. - The startDate or endDate is in the future.	HTTP 400	1002001
- The client is not authorized to make this request. - The requested patient was not accessible or was not found.	HTTP 401	1002001

Example:

```
{
  "Response": {
    "responseCode": "1002001"
  }
}
```

Note: Data is fictional.



ResMed Corp

9001 Spectrum Center Boulevard San Diego CA 92123 USA

See www.resmed.com for other ResMed locations worldwide.

For patent and other intellectual property information, see www.resmed.com/ip.

AirView is a trademark and/or registered trademark of the ResMed family of companies. 2281511/3 2021-06

ResMed.com